**EXPANDABLE
CORPECTOMY
CAGE**

**EXPANDABLE
CORPECTOMY CAGE**
GENİŞLEYEBİLİR KORPEKTOMİ KAFES

The key in

orthopaedic

sciences

**SURGICAL TECHNIQUE
CERRAHİ TEKNİK
IMPLANTS & INSTRUMENTS
İMLANTLAR VE EL ALETLERİ**

CONTENTS İÇİNDEKİLER

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Introduction

Giriş

EN The Expandable Corpectomy Cages may be used after complete or partial destructions of the vertebral body. Such destructions may be caused by;

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EN Expandable Corpectomy Cage is applied if the anterior thoracic or lumbar spine needs to be reconstructed.

Features of the Expandable Corpectomy Cage;

- One or more vertebral bodies can be substituted
- Expandable Corpectomy Cage may be distracted continuously in situ
- Expandable Corpectomy Cage corrects malalignments of the spine by distraction
- Expandable Corpectomy Cage is anchored easily and safely between the adjacent end plates of the vertebral bodies.
- A small number of instruments makes handling easy.
- The Corpectomy Cage has a hollow design. If necessary the implant may be filled with bone graft through large orifices in the cylindrical portions. This may be done either before or after insertion.

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EN The Expandable Corpectomy Cage consists of three components: a center piece and two end pieces. The implant is available in three external diameters, 20 mm and 24 mm. The center piece's one side is equipped with a right and the other side is equipped with a left thread. Accordingly, the end pieces are equipped with either a left or a right thread which are marked with "R" and "L".

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EN The maximum and minimal distraction length depend on the center piece. To avoid overdistracton the section coloured in blue indicates the maximum distraction height . The distraction heights range from 25-34 mm to 74-114 mm. Spikes at the rim of the end pieces improve the anchorage of the implant on the end plates of the vertebral bodies.

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EN The Unilateral Expandable Corpectomy Cage has already been used for some years as a vertebral body replacement for the anterior cervical and thoracic vertebral column. Unilateral Expandable Corpectomy Cage can also be used for intersomatic fusion.

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EN Unilateral Expandable Corpectomy Cage is used for the reconstruction of substance defects, support of the anterior column as well as restoration of the vertebral column profile. An additional anterior or posterior stabilization by means of a rod/screw system or a plate system is necessary.

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EN Unilateral Expandable Corpectomy Cage is available in four outer diameters (10, 12, 14 and 16). The distraction ranges extend from 7-10 mm to 39-65 mm. In part, the Unilateral Expandable Corpectomy Cages are available with an angle of 6°.

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EN By turning the distraction ring the implant is spread open and the distraction height becomes fixed. The distraction height of the Unilateral Corpectomy may be adjusted continuously as per the individual case.

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**Material
Malzeme**

EN Certified Ti6Al4Veli alloy according to ASTM F 136 - ISO 5832-3. It is selected for its, modulus of elasticity closest to that living bone, high fatigue strength and perfect biocompatibility.

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Please Note : This document is intended as a guide for the surgeon only. There are multiple techniques for the insertion of Expandable Corpectomy Cage and, as with any surgical procedure, a surgeon should be thoroughly trained and beware that this procedure is appropriate for the patient before proceeding.

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Exposure Açma İşlemi

EN Determine the correct approach. The involved segment(s) of the vertebral column are exposed. Complete or partial resection of the vertebral body (bodies) including the adjacent discs is performed.

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EN The gap to be bridged is measured using the measuring caliper (Ref:1068306001). Corresponding size with measuring calliper on the measuring template (Ref:1083070001) shows the appropriate size of the implant.

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⚠ Attention / Dikkat

EN The first section on the measuring template corresponds to the Expandable Corpectomy Cage that has not been distracted. The second(striped) section illustrates the possible distraction height.

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
EN Unilateral Expandable Corpectomy Cage's measuring table is below the Expandable Corpectomy Cage's table on the measuring template. Chose the appropriate Unilateral Expandable Corpectomy Cage to the measuring caliper.

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

Distraction Distraksiyon

EN The expansion instrument (Ref:10683010005) is inserted into the rotation holes of the center piece. Distraction is achieved by turning the center piece as indicated by the arrows.

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! Attention / Dikkat

EN As soon as the blue portion becomes visible, the maximum distraction height is achieved.

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EN As an alternative to the Expansion Instrument the Expansion Instrument/5 may be used. The expansion instrument/5 is attached to the rotation hole of the center piece of the cage. To achieve desired distraction, expansion instrument/5 turned in the direction of the arrow.

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EN As soon as the appropriate distraction height has been achieved the fixation instrument (Ref:10683050001) is removed.

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**Attention / Dikkat****EN**

Check fitting and positioning of the implant.

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
**Locking the Implant
İmplantın Kilitlenmesi****EN**

After the correct position of the implant has been checked the end pieces are locked to the center piece using a set screw (Ref:10133202010). Again the screw driver for set screw (Ref:10683030001) is used.

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**Bone Grafting
Kemik Grefti****EN**


The void in the Expandable Corpectomy Cage resulting from the distraction can be topped up through the large drill holes in the end pieces. Additional bone grafting around the implant is possible. The Expandable Corpectomy Cage may also be filled with bone cement depending on the situation.

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

Unilateral Expandable Corpectomy Cage

10132102100	Ø10 - 7 X 10			10232122606	Ø12 - 39 X 65	6°
10232102200	Ø10 - 10 X 13			10232142100	Ø14 - 10 X 13	
10232102300	Ø10 - 12 X 17			10232142200	Ø14 - 12 X 17	
10232102400	Ø10 - 16 X 25			10232142300	Ø14 - 16 X 25	
10232102406	Ø10 - 16 X 25	6°		10232142306	Ø14 - 16 X 25	6°
10232102500	Ø10 - 24 X 40			10232142400	Ø14 - 24 X 40	
10232102506	Ø10 - 24 X 40	6°		10232142406	Ø14 - 24 X 40	6°
10232102600	Ø10 - 39 X 65			10232142500	Ø14 - 39 X 65	
10232102606	Ø10 - 39 X 65	6°		10232142506	Ø14 - 39 X 65	6°
10232122200	Ø12 - 10 X 13			10232162100	Ø16 - 10 X 13	
10232122300	Ø12 - 12 X 17			10232162200	Ø16 - 12 X 17	
10232122400	Ø12 - 16 X 25			10232162300	Ø16 - 16 X 25	
10232122406	Ø12 - 16 X 25	6°		10232162306	Ø16 - 16 X 25	6°
10232122500	Ø12 - 24 X 40			10232162400	Ø16 - 24 X 40	
10232122506	Ø12 - 24 X 40	6°		10232162406	Ø16 - 24 X 40	6°
10232122600	Ø12 - 39 X 65			10232162500	Ø16 - 39 X 65	

Instruments / El Aletleri

10683010005		Expansion Instrument
10683020002		Expansion Instrument / 2
10683020003		Expansion Instrument / 3
10683020005		Expansion Instrument / 5
10683030001		Screw Driver for Setscrew
10683040001		Impactor / Straight
10683040002		Impactor / Angled
10683050001		Fixation Instrument
10683060001		Measuring Caliper
10683070001		Measuring Table

The products being manufactured by TIPSAN A.S. has been certified from SZUTEST (Strojirensk•
Zkušební Ustav) for EC certificate **(CE 1015)** according to Annex II. 3 and II.4 of Council Directive
93/42/EEC concerning medical devices

Kemalpasa Cd. 7404/1 Sk. No: 3 Pinarbasi - IZMIR - TURKEY
Tel: +90 232 479 56 54 Fax: +90 232 479 58 27
marketing@tipsan.com.tr www.tipsan.com.tr