
**MODULAR STRAIGHT STEM
C E M E N T E D**

**MODULAR STRAIGHT
STEM CEMENTED**

The key in

orthopaedic

sciences

**SURGICAL TECHNIQUE
CERRAHİ TEKNİK
IMPLANTS & INSTRUMENTS
İMLANTLAR VE EL ALETLERİ**

CONTENTS İÇİNDEKİLER

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Introduction

Giriş

EN The modular straight stem hip prosthesis is designed to be used as a primary and revision total hip. The 4 stem sizes allow accurate matching of prosthesis to the medullary canal. The prosthesis offers dual fixation by combining the locking effect of a stem to bone contact, plus minimal bone cement support.

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
Material Malzeme

EN Certified Stainless Steel according to ASTM F 138 and ISO 5832-1. Available in CoCr alloy according to ASTM F 75 and ISO 5832-4.

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Please Note : This document is intended as a guide for the surgeon only. There are multiple techniques for the insertion of Modular Straight Stem Cemented and, as with any surgical procedure, a surgeon should be thoroughly trained and beware that this procedure is appropriate for the patient before proceeding.

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Femoral Neck Osteotomy Femoral Boynun Osteotomisi

EN Because of the normal anteversion of the femoral neck, the prosthesis can only be positioned in the correct axial direction if the length of neck stump is medially no more than 1 cm.

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EN A very important step is the trapezoidal opening of the medullary cavity. This opening, which may be made with a Tuke saw, should reach the tip of the greater trochanter. Then the femoral head is removed using a Femoral Head Extractor (Ref: 10606100001). The resected cancellous bone is removed with a gouge to expose the medullary canal.

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EN The awl (Ref:10610331003) is introduced in the region of the resected tip of the greater trochanter and should follow the long axis of the femur. It will orientate itself to the course and shape of the medullary canal and provide a guideline for the ideal seating of the stem.

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
EN Remove any loose, unsupportive, cancellous bone from the canal with a Curettes (Ref: 10605210005-10605210025-10605210103). Close the distal canal with a medullary plug at least at 1 cm distal to the tip of the stem. Pay attention in the choice of the distal plug, to ensure that it can resist to the cement pressure. Clean the intramedullary canal with pulse lavage and dry it. Keep the canal packed until cement is ready to be injected.

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EN Introduce the cement into the canal retrograde using a cement gun. Pressurize the cement column to allow the cement to interdigitate into the cancellous bone.

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EN The stem is introduced using the impactor (Ref: 10607031321), until prosthesis's colar will reach the same level with the resection level.

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EN After the cement has hardened completely, and the taper has been carefully rinsed, the ball head is mounted.

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EN The definitive head is positioned with a slight rotational movement. Seat the head tightly with tap on the Femoral Head Impactor (Ref:10607020003).

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EN Reduction and check of functioning.

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Attention/Dikkat

EN

Use the Impactor-Extractor Pim (Ref: 10607031321) and extracting hole of the prosthesis for removing the prosthesis

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


The products being manufactured by TIPSAN A.S. has been certified from SZUTEST (Strojirensk• Zkušební Ustav) for EC certificate **(CE 1015)** according to Annex II. 3 and II.4 of Council Directive 93/42/EEC concerning medical devices

Kemalpasa Cd. 7404/1 Sk. No: 3 Pinarbasi - IZMIR - TURKEY
Tel: +90 232 479 56 54 Fax: +90 232 479 58 27
marketing@tipsan.com.tr www.tipsan.com.tr