PRESS FIT CUP UNCEMENTED

NEW

PRESS FIT CUP

The key in

orthopaedic

sciences

**SURGICAL TECHNIQUE
CERRAHI TEKNİK
IMPLANTS & INSTRUMENTS
İMLANTLAR VE EL ALETLERİ**

CONTENTS İÇİNDEKİLER

Introduction

Giriş

2

Reaming

Oyma

5

Cup Trialing

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6

Impacting of the Acetabular Shell

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8

Screw Fixation

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9

Insert Impaction

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10

Instruments&Trays

El Aletleri&Tavalar

12

Introduction

Giriş

Press Fit Cup

Pres Fit Kap

EN In principal for all forms of arthrosis, advance wear and tear of the hip joint because of degeneratif, posttraumatic or rheumatoid arthritis fracture of avascular necrosis of the head of femur. Sequelae of earlier operations such as osteosynthesis, joint reconstruction, arthrodesis, hemiarthroplastic or total hip replacement. Press-fit cup, is an acetabular implant for uncemented fixation. The titanium cup has three holes on supero-lateral side. The screws have a spheric head and have always a good contact with their seatings, even if the introduction angle may be different. This prevents the fracting corrosion and assures the screw stability in spite of the different possible introduction angles of the screw.

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Material

Malzeme

EN Certified Titanium alloy according to ASTM F 136 - ISO 5832-3. It is selected for its, modulus of elasticity closest to that living bone, high fatigue strength and perfect biocompatibility.

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Press Fit Cup Insert Pres Fit Kap İnsert

EN The press fit cup can reduce the possibility of wear debris from the inner surface of the UHMWPE cup. It is done by machining and polishing of the inner surface. The good precision of manufacture assures more over the direct contact between the UHMWPE insert and inner surface of the metal cup reducing micromotions and fretting.

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Material Malzeme

EN PE-UHMW crosslinked material meets ASTM F 648 and ISO 5834-1/2 standards. Selected for its, high degree of purity, good biotolerance, good mechanical performance and friction properties.

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Starex Screw

Starex Vida

EN If bone quality is such that the titanium cup has to be secured with screws, special bone screws with spheric head and cancellous thread (\varnothing 6.5 mm) are used for fixation. For fixation, flexible drill and universal joint screw driver are available. The screw have always a good contact with their seatings, even if the introduction angle may be different. This prevent the fretting corrosion and assure the screw stability.

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
Material

Malzeme

EN Certified Titanium alloy according to ASTM F.136 and ISO 5832/3.

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
Please Note : This document is intended as a guide for the surgeon only. There are multiple techniques for the insertion of Press Fit Cup Uncemented and, as with any surgical procedure, a surgeon should be thoroughly trained and beware that this procedure is appropriate for the patient before proceeding.

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Reaming Oyma

EN After arthrotomy and osteotomy of the femoral neck, expose and prepare the acetabular cavity, remove osteophytes.

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
EN Start reaming with acetabular reamers. The ideal reaming axis has an inclination of 40°/45°, and an anteversion of 15°/20°. (Anteversion recommended for posterior approaches).

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EN Reaming of the acetabulum starts with the smallest reamer and increases in increments of 2 mm, until a perfectly regular hemispherical cavity has been obtained.

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Attention/Dikkat

EN During final reaming, avoid changing the reamer axis, in order not to make the reaming oval, which may affect or prevent primary seating.

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
EN As a general rule the right diameter corresponds to 4 or 6 mm greater than the femoral head diameter size. Take care to retain, as far as possible, the bone stock up to the level of anterior and posterior columns. Reamings may be saved for void filling between implant and acetabulum.

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EN Assemble the trial cup with the same diameter of the last reamer with the handle. Insert the trial cup into the reamed cavity in order to estimate the depth and the orientation of the acetabular component.

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EN Trial cups:

- Are smooth and have the same dimensions as the reamers to avoid damaging the socket.
- Are slightly undersized compared to the implant to allow a maximum press-fit effect with the definitive implant.
- Have several openings to permit a direct view of the underlying acetabular surface.

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EN Both the implant and trial cup have a 5° raise. Marks on the trial cup or acetabular shell help to identify coverage top (see image).

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i Additional info/Ek bilgi


EN In order to ensure the definitive acetabular shell correct positioning in depth, mark the coverage top. A positioning guide is available to aid in the positioning.

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! Attention/Dikkat

EN If the trial cup is not stable or primary stability is doubtful, especially in the presence of poor bone quality, it is possible to choose a larger cup size, either with or without additional acetabular reaming.

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EN After a satisfactory trial, acetabular shell can be positioned. The definitive acetabular shell size will be the same as the final reamer size.

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
i Additional info/Ek bilgi

EN A positioning guide is available to aid in the acetabular shell positioning: the positioning guide will be positioned on the top of the press-fit cup handle. The two rods of positioning guide are inclined at 45° and 20° to the press-fit cup handle.

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
EN The assembly is positioned in the acetabular shell.

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EN Impact the acetabular shell with the aid of a hammer, until fully seated and stabilized. Note that the final acetabular shell is press fit, not threaded into the acetabulum.

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EN Should additional screw fixation be necessary, holes are drilled (Ref:10809071040-50) in the direction of the sacro-iliac joint using the drill guide (Ref:10608021021) and a drill with flexible drive shaft (Ref:10601361101).

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EN The depth of the drilled hole in the bone is measured using the screw depth gauge (Ref:10608031021), in order to determine the correct screw length. Following this, the first thread pitches (about 10mm deep) for the spongiosa screw are pre-cut using the screw tap (Ref:10606001101) on a Universal Joint Tap Handle. The screw is then easily inserted into the deeper soft spongiosa layers.

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EN Where high torque is necessary for insertion, it is recommended to either revert to a shorter screw or to pre-cut the thread in the bone to the measured length using a steel screw. Following this, the steel screw is removed and replaced with a titanium screw.

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EN When the screw holes are adequately prepared, the screws are inserted in the direction of the sacro-iliac joint and tightened using the Universal Joint Screw Driver (Ref:10604010012).

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EN Clean the interior surface of the acetabular shell. Position the trial liner corresponding to the diameter foreseen. If a PE Insert will be used, stability tests are performed after having positioned the trial or final stem and the trial head. Reduce the hip in order to test the joint stability and limb length.

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! Attention/Dikkat

EN Test of stability must be performed with trial heads and not with final heads. The head sizes XL and XXL (for Ø28 mm) have a skirt. This may decrease the Range of Motion and may cause an impingement risk with the liner.

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
EN When the screw position is correct, the flat screw heads are fully recessed into the screw sockets, avoiding contact with the polyethylene insert.

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
EN The polyethylene insert has two small holes on its front face, which grip the pins of the Screw-in Attachment (Ref:10601361028).

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
EN First, the polyethylene insert is attached to impactor for PE liner (Ref:10601361200-10-15-20) then assembled to impactor handle. Insert placed into titanium shell with aligning the lines and fixed with pressure. The raised rim of the insert, which serves as protection against dislocation in extreme positions without load especially in low sitting positions with adducted thigh, is normally inserted in a cranio-dorsal position, i.e. at an approximately "10 o'clock" position for the right hip, and "2 o'clock" for the left hip. When the desired position is achieved, the polyethylene insert is impacted into the titanium shell using the hammer.

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

Ø3.5x30x50

Notes/Notlar

The products being manufactured by TIPSAN A.S. has been certified from SZUTEST (Strojirensk•
Zkušební Ústav) for EC certificate **(CE 1015)** according to Annex II. 3 and II.4 of Council Directive
93/42/EEC concerning medical devices

Kemalpasa Cd. 7404/1 Sk. No: 3 Pinarbasi - IZMIR - TURKEY
Tel: +90 232 479 56 54 Fax: +90 232 479 58 27
marketing@tipsan.com.tr www.tipsan.com.tr