

PFSN
PROXIMAL FEMORAL
SHORT NAIL

PFSN
PROXIMAL FEMORAL SHORT NAIL
PROKSİMAL FEMUR KISA ÇİVİ

The key in

orthopaedic

sciences

SURGICAL TECHNIQUE
CERRAHI TEKNİK
IMPLANTS & INSTRUMENTS
İMLANTLAR VE EL ALETLERİ

CONTENTS İÇİNDEKİLER

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Please Note : This document is intended as a guide for the surgeon only. There are multiple techniques for the insertion of Proximal Femoral Short Nail System and, as with any surgical procedure, a surgeon should be adequately trained and be aware that this procedure is appropriate for the patient before proceeding with this technique.

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-Distal diameters; Ø10, Ø11, Ø12, Ø13
-Length; 180mm, 200mm, 220mm

-Distal çap; Ø10, Ø11, Ø12, Ø13
-Uzunluk; 180mm, 200mm, 220mm

Locking Screws - Captured Ø5mm
Distal Kilitleme Vidaları, Ø5mm

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
EN Position patient supine on an extension table or a radiolucent operating table. Position the C-arm of the image intensifier in such a way that it can visualize the proximal femur exactly in the lateral and AP planes.

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EN For unimpeded access to the medullary cavity, abduct the upper part of the body by about 10–15° to the contralateral side (or adduct the affected leg by 10–15°).

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EN Under image intensifier control, place the Ruler (Ref: 10682100003) on the femur and position the square marking over the isthmus. If the transition to the cortex is still visible to the left and right of the marking, the corresponding nail diameter may be used. If the PFSN Ø 10.0 mm is still too thick, the medullary canal must be reamed accordingly.

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**Attention/ Dikkat**

EN Excessive force will have to be used if too thick a nail is introduced. This may result in loss of reduction and bone fractures.

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EN In the AP view, the nail insertion point is normally found on the tip or slightly lateral to the tip of the greater trochanter in the curved extension of the medullary cavity.

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
EN The mediolateral angle of the implant is 5°. This means that the 2.5 mm Guide (Kirschner) Wire must be inserted laterally at an angle of 5° to the shaft. On the lateral view, the guide wire must be located in the centre of the medullary canal, up to a depth of 15 cm. The guide wire can be inserted either manually with the Universal Chuck with T-Handle or a power tool with the Quick Coupling for Kirschner Wires.

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⚠ Attention/ Dikkat

EN A correctly selected insertion point and angle are essential for a good surgical result. Check the correct position of the Kirschner (guide) wire with the aid of x-rays and by positioning a nail anteriorly on the femur.

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EN Guide the Cannulated Drill (Ref: 10682050017) through the protection sleeve (Ref: 10682061722) over the guide wire and ream manually with the Universal Chuck with T-Handle as far as the stop on the protection sleeve.

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⚠ Attention/ Dikkat


EN One should be very careful while drilling unstable, multifragmentary fractures. You should avoid varus displacement of the medial fragment by making sure that the hole is drilled both in the medial fragment and the lateral part of the femur.

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Option: Entry with Kuntscher BIZI Alternatif Yol: Kuntscher BIZI ile giriş

EN Open the femur or enlarge the entry point with the Kuntscher BIZI. Drive the Kuntscher BIZI over the guide wire into the femur until the curved part on the Kuntscher BIZI shaft is level with the trochanter tip.

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
Assemble Instruments El Aletlerinin Montajı

EN Guide the Guide Screw (Ref: 10682090004) through the Nail Holder (Ref: 10682090001) and secure the nail tightly to the Nail Holder using the Screw Driver / 8 (Ref: 10682010008). The nail diameter has already been determined during preparations for surgery. Ensure that the connection is tight to avoid deviations when inserting the screws through the guide (Ref: 10682090003). Do not attach the guide yet.

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EN Carefully insert the nail manually as far as possible into the femoral opening. Slight twisting hand movements may help the insertion. If the nail cannot be inserted, select a smaller nail diameter.

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EN Insertion can be supported by light blows with the synthetic Hammer on the mounted Nail Holder Attachment (Ref: 10682090002) of the Nail Holder (Ref: 10682090001).

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EN The correct nail insertion depth is reached if the future position of the proximal locking screw is just above the calcar in the distal half of the femoral neck. The future position of the proximal locking screw can be assessed on the AP view by means of the position of the proximal holes in the nail. Ensure that sufficient space is also available for the other proximal screw.

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Attention / Dikkat

EN – If the medullary canal is too narrow, it should be reamed to at least \varnothing 11 mm.
– Make sure that the nail is tightly screwed to the Nail Holder before inserting the nail.

– If a hammer to be used, tap gently on the Nail Holder. The use of excessive force will result in loss of reduction or secondary fractures. Apply taps only to the Nail Holder Attachment.


– An excessively cranial or caudal nail position will result in incorrect placement of the screws and should therefore be avoided at all costs.

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EN Secure the corresponding Guide tightly to the Nail Holder (Ref: 10682090001). Select the system consisting of protection sleeve (Ref: 10682120008), drill sleeve (Ref: 10682072509) and trocar (Ref: 10682030190).

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i Additional Info / Ek Bilgi

EN The position of the nail can now be checked on the AP view by positioning a guide wire over the Guide.

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EN Make a incision and insert the selected Drill Sleeve System through the guide until the bone is reached. Mark the femur and remove the trocar.

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
EN Insert a new 2.5 mm Threaded Kirschner Wire (Ref: 10682110001) through the proksimal drill sleeve, check direction and position under the image intensifier in AP and lateral views. The correct position of the proksimal locking screw tip on the AP image is at a distance of 5 to 10 mm from the subchondral bone. The Kirschner wire will therefore need to be inserted to the subchondral bone or up to a maximum distance of 5 mm away. In the lateral view the wire should be located in the centre of the femoral neck.

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! Attention / Dikkat

EN If the Kirschner (guide) wire is not in the desired position or if it is bent it must be reinserted. If the nail positioned incorrectly, remove the Kirschner wire and the drill sleeve system. The nail must be repositioned by means of rotation, deeper insertion or slight retraction. Reset the drill sleeve system and introduce a new guide wire.

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
EN Make sure the Nail Holder is not moved before inserting the lower proximal locking screw guide wire. Drilling over a bent wire guide can result in drill breakage or damaging of the nail itself.

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EN Insert Drill Sleeve System through the upper proximal locking hole on the guide to the bone. Then remove the trocar and insert a second, new 2.5 mm Kirschner wire (Ref: 10682110001) through the drill sleeve (Ref: 10682072509) into the bone. The insertion depth of the Kirschner (guide) wire should be 10 mm less than the insertion depth of the lower proximal locking screw Kirschner (guide) wire. This ensures that the lower proximal locking screw will not take weight.

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! Attention / Dikkat


EN Verify that the kirschner wires are parallel in both planes and that their tips form a horizontal line on the AP view.

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Remove the Drill Sleeve Drill Sleeve'in Sökülmesi

EN Drill sleeve (Ref: 10682072509) must be carefully removed without moving the Kirschner wires from their original positions. Before the length measurement, recheck the position of the Kirschner wires on the AP view.

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
EN In order to prevent possible rotation of the medial fragment when inserting the upper proximal locking screw, prior insertion of the lower proximal locking screw is recommended.

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EN Guide the Measuring Device (Ref: 10682100001) into the protection sleeve (Ref: 10682060911) to the bone and determine the required length of the proximal screw. The length of this screw is indicated on the measuring device and is calculated to end 5 mm before the tip of the guide wire.

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
EN Assemble the Drill Stopper (Ref: 10682050001) to Reamer for proksimal screw (Ref: 10682060008) with considering length of measured nail and fix it.

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EN Carry on drilling till drill stopper (Ref: 10682050001) is stopped by protection sleeve (Ref: 10682060911). As the tip of the proksimal locking screw is selftapping, usually no further drilling and tapping is needed.

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Procedure with Hard Bone Sert Kemikteki İşlemler


EN With hard bone, further drilling and tapping with the Tap is recommended up to the length of the upper proksimal locking screw previously measured.

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EN Assemble the upper proksimal locking screw to the screwdriver for proksimal screw, insert the selected upper proksimal locking screw over the guide wire completely to the stop. Remove and discard the guide wire (Ref: 10682110001) of the upper proksimal locking screw.

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
! Attention / Dikkat

EN Do not insert the upper proximal locking screw with undue force. Ensure that the lateral end of the upper proximal locking screw clearly protrudes from the lateral cortex and is not inserted into it.

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
EN To Measure length of the Lower proximal screw, to drill hole and to insert screw operations are performed in accordance with operations performed for upper proximal screw.

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
Distal Locking
Distal Kilitleme

EN Choose static or dynamic distal locking. Select the appropriate targeting module hole that corresponds to the nail size implanted. Locking options include static, dynamic, both or none.

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EN Distal locking is usually performed with a single locking bolt. For static interlocking use the caudal locking hole only; for dynamic interlocking the cranial locking hole. Subtrochanteric fractures may be double-locked. Secondary dynamization is possible by postoperative removal of the static locking bolt.

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EN Make a incision and insert the green Drill Sleeve System [protection sleeve (Ref: 10682060911)+drill sleeve distal (Ref: 10682074009)+trocar (Ref: 10682030185)]through the locking hole selected in the guide until the bone is reached.

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EN Remove the Trocar and drill through both cortices using the Drill. Determine the length of the required locking with screw length gauge. Add 2 to 4 mm to the reading to ensure that the thread engages the far cortex. If requested screw length can be determined by reading the scale on the drill during drilling operation.

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! Attention / Dikkat

EN – Before distal locking, ensure that no diastasis has occurred intraoperatively. Healing may be delayed if the distal locking bolt are secured despite the presence of diastasis.
– Ensure that all connections between the nail, PFSN Holder and guide are still secure, otherwise the nail itself may be damaged during drilling of the distal locking holes.

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
EN Insert the locking screw through the protection sleeve using the screwdriver for distal screw (Ref: 10682010005) and finish distal locking.

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
EN Align the end cap with the nail axis using the screw driver for distal screw (Rev:10682010005) in order to avoid unsuitable assembly of end cap. Screw the end cap completely onto the nail until it stops.

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

Screw Driver for Proximal Screw
Proksimal Vida Tornavidası

10682010005

Screw Driver for Distal Screw
Distal Vida Tornavidası

10604051003


Double Open Wrench
İki Ağızlı Anahtar

10682020001

Compression Device for Proximal Screw

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Drill Sleeve - Proximal
Proksimal Dril Sleeve

10682074009

Drill Sleeve - Distal
Distal Dril Sleeve

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10682090005

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10682110001

Threaded Kirschner Wire
Ucu Yivli Kirschner Teli

10682050004

Dril - Calibrated
Kalibreli Drill

The products being manufactured by TIPSAN A.S. has been certified from SZUTEST (Strojirensk•
Zkušební Ústav) for EC certificate **(CE 1015)** according to Annex II. 3 and II.4 of Council Directive
93/42/EEC concerning medical devices

Kemalpasa Cd. 7404/1 Sk. No: 3 Pinarbasi - IZMIR - TURKEY
Tel: +90 232 479 56 54 Fax: +90 232 479 58 27
marketing@tipsan.com.tr www.tipsan.com.tr